



AMMONIUM TEST KIT FOR WATER

(Sensitivity 0,5 to 10 ppm NH₄₎

Product Code: SP-KS-

- 1. Rinse the test tube with the water to be tested and fill to the 2,5 ml mark and complete to 5 ml mark with **Ammonium-A Indikator**, place stopper and shake well 10 minutes.
- 2. Add 3 drops of **Ammonium-B Indicator**, then add 3 drops of **Ammonium-C Indicator**.
- 3. Place stopper and shake well.
- 4. Compare the color with the color reference chart. Same color indicates ammonium in ppm(mg/l) units.

Information:N=NH₄ x 0,77